

MedReady, Inc
2103 Border Avenue
Torrance, CA 90501

Phone: 310.328.7557
Fax: 310.328.7773
e-mail: sales@medready.net
Web: www.medreadyinc.net



Thank you for choosing MedReady Automated Medication Dispenser.

Please read through the operating instructions so you will know how to operate your MedReady properly. If you have any questions regarding the operation of your unit, please call MedReady Tech Support at the number above, or send us an email.

CONTENTS

1 Medication Dispenser (incl. 1 tray) 1 Rechargeable NiMH Battery Pack 2 Keys
2 Emergency Information Cards ** 2 Loading Rings ** 1 set Instructions
1 AC adapter (120v, 60 hertz) 1 Telephone cable connector (1750 modem only) **

** - Not included with 1700PRN

NOTE: For convenience, you can store the Emergency Cards and the Loading Rings under the clip beneath the device

IMPORTANT NOTES

1. Only pills should go in the tray. NITROGLYCERINE, liquids and powders **SHOULD NOT BE PLACED IN THE DEVICE**. As with all medication, please keep away from children. The tray should be cleaned regularly.
2. MedReady is powered by **BOTH** its AC adapter and its battery pack. Please plug in the battery pack **FIRST** and then plug in the AC adapter. The AC adapter must charge the battery pack for at least 4 hours before MedReady's first use.
3. You should use the AC wall power whenever possible. The back-up batteries should only be used without A/C power for short-term power outages and short-term travel of less than 72 hours. When the device is using battery power, **ONLY** the two dots (:) will be displayed by the clock.
4. The key **CANNOT** be removed when the lock is in the open position. This is to remind you that the MedReady is not locked and also to help you keep track of the key.
5. For maintenance purposes, you should replace the battery pack once a year. Use only MedReady battery packs.

PLEASE READ AND UNDERSTAND THE FOLLOWING:

MedReady is designed as a medication aid only. MedReady is NOT a medication Safe or Vault, and it is NOT TAMPERPROOF. With determination, a patient will be able to break into the device and access more medication than intended. If there is any doubt a patient will have an uncontrollable urge to access more medication than prescribed, please do NOT use this device.



LIMITED WARRANTY

MedReady, Inc. warrants this product is free from original manufacturing defects in material or workmanship under normal use for one (1) year from the date of purchase. It does not cover damage resulting from accident, misuse, abuse, neglect, mishandling misapplication or alteration of this product or its documentation. MedReady's liability is limited to the cost of repair or replacement, at our option, of any defective part or product and shall in no event include incidental or consequential damages of any kind.

If the product proves defective during this warranty period, MedReady will pay all labor and material expenses for the covered item but the customer must pay any shipping charges to return the item to MedReady. If the repair is covered by the warranty, MedReady will pay the return shipping charges.

MedReady reserves the right to revise or update its products or documentation without obligation to notify any individual or entity. MedReady makes no warranty or representation, expressed, implied or statutory, with respect to its products, or the contents or use of this documentation, and specifically disclaims its quality, performance, merchantability, or fitness for any particular purpose.

MedReady Replacement parts can be ordered at www.medreadyinc.net/

2302	Blue top with door without lock	\$12.00
2302-L	Blue top with door and lock and 2 keys	18.95
5300	Replacement gear assembly	9.00
5310	Tray and clear plastic cover for pre-filling	18.00
5313-12	Replacement tray	9.00
5313-13	Replacement clear plastic cover	9.00
6308-K	Replacement Key #1688	4.00
6311-00	Replacement battery	12.85
6312	Replacement AC adapter	8.50
5307-12	Battery Door	1.20

Units can be refurbished for approximately \$70.00 including standard return shipping.