

# MedReady™

Available with cellular  
compliance monitoring

- *Timed-alarm system provides next dosage automatically*
- *Holds up to 7 days of 4 times a day dosing*
- *MedReady's functions and capabilities provide the best value in the medication dispensing marketplace*
- *The cellular dispenser has the capability to communicate compliance events via cell towers. No phone line or internet connection is necessary at the dispenser location.*
- *Check compliance history, change dispensing times, and update all device settings via the [www.medready.net](http://www.medready.net) website.*



## Automated Medication Dispenser

[medreadyinc.net](http://medreadyinc.net)

MedReady, Inc.  
2103 Border Ave.  
Torrance, CA 90501  
310-328-7557  
[Sales@MedReady.net](mailto:Sales@MedReady.net)

*One Less Argument • More Independence  
Feel Better • Great Price*





## FEATURES



Large 1/2" L.E.D. displays the time and aids setting alarms.

- Medication regimen can dose up to 4 times per day:
  - ◆ "Time for medicine" alarm beeps 5 seconds on, 1 second off continuously until the medication door is slid open and the "off" position is activated.
  - ◆ Note: Disabling alarm requires the medication door to be slid open.
  - ◆ Programmable alarm sounds for 3, 10, 20, 30, ...250 minutes, as user selects, and can chirp each minute thereafter until the dose is taken or the next one is available.
  - ◆ Powerful, 90dB, 400Hz. audible alarm is designed for maximum effectiveness.
  - ◆ Audible and visual alarms can be independently turned on or off.
- Only the most recent medication dose is available. Pill tray rotates to the next dose when the next audible alarm sounds. If a patient misses a dose, the pill dispenser denies access to the missed dose.

# MedReady™

Automated Medication Dispenser  
Available with cellular  
compliance monitoring

- MedReady holds 28 doses of dry medication such as pills or capsules.
  - ◆ 4 Doses a day for 7 days
  - ◆ 3 Doses a day for 9 days
  - ◆ 2 Doses a day for 14 days
  - ◆ 1 Dose a day for 28 days
- Alarm in pill dispenser is as easy to set as a digital clock.
- Medication is under lock and key, and is tamper resistant.
- Pill compartments are "ramped" for easy removal and can fit up to 9 regular **m&m's**.
- Compartment divisions keep medications in place even if the unit is turned upside-down.
- MedReady should remain plugged into electric power, although there is a 48 hour rechargeable battery back-up for travel or power failures. MedReady should not be used in wet or very damp environments.
- "OFF" feature for long term disabling of the dispenser.
- "Early dosage" feature allows removal of medication prior to scheduled times.

- Cellular Communication from the MedReady dispensers to our servers.
  - ◆ Compliance notifications via email, text, and automated phone calls.
  - ◆ All medication compliance events are communicated in real time.
  - ◆ Monthly service plan includes monitoring and MedReady.net compliance website access.
- Manage MedReady online
  - ◆ Set dispensing times, change device settings, trigger early doses view detailed reports.

Remaining Doses: 24

Days Remaining: 6

ALARM TIME: Friday 11/06/15

ALARM TIME	DOSE TIME
06:00	06:10
09:00	09:05
15:00	16:10
18:00	18:15

Alarm	Time	Active?	Meds	Action
⊙ Alarm 1	6:00 AM	Active: <span style="color: green;">Yes</span>		<a href="#">Edit</a>
⊙ Alarm 2	9:00 AM	Active: <span style="color: green;">Yes</span>		<a href="#">Edit</a>
⊙ Alarm 3	3:00 PM	Active: <span style="color: green;">Yes</span>		<a href="#">Edit</a>
⊙ Alarm 4	6:00 PM	Active: <span style="color: green;">Yes</span>		<a href="#">Edit</a>

Other Functions

Alarm Test Send Early Dose Send Cancel Early Dose Send

## OPTIONS

- 1. MedReady Monitoring (MR-357)** The user does not need a cell phone, MedReady uses its own Internal Cellular Device
- 2. MedReady Monitored1750** using Plain Old Telephone Service (POTS)
- 3. FL - Alarm With Flashing Light** for hearing impaired. Audible alarm PLUS Flashing Light.
- 4. 5310 Extra Pill Tray** with clear plastic lid. Preload an additional tray for convenience.

